

IHE Integration Statement			
Vendor	Product Name	Version	Date
InsiteOne LLC	Unified Clinical Archive	11.1	8/20/2022

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profile	Actor	Integration Profile Option
Access to Radiology Information	Image Manager/Archive	None
Access to Radiology Information	Report Repository	None
Cardiac Catheterization Workflow	Image Manager/Archive	None
Cardiac Catheterization Workflow	Image Manager/Archive	Cardiac Cath
Cardiac Catheterization Workflow	Image Manager/Archive	Intermittently Connected Modality
Cardiac Catheterization Workflow	Performed Procedure Step Manager	None
Consistent Presentation of Images	Image Manager/Archive	None
Consistent Time	Time Client	None
Digital Breast Tomosynthesis	Image Manager/Archive	None
Digital Pathology Image Acquisition	Image Manager/Archive	None
Echocardiography Workflow	Image Manager/Archive	Echocardiography
Echocardiography Workflow	Image Manager/Archive	Intermittently Connected Modality
Echocardiography Workflow	Performed Procedure Step Manager	None
Evidence Documents	Image Manager/Archive	None
Cardiology Evidence Documents	Image Manager/Archive	None
Image Object Change Mgmt	Image Manager/Archive	None
Key Image Note	Image Manager/Archive	None

Internet address for vendor's IHE information:

https://insiteone.com/s/UCA_R11_1_IHE_Integration_Statement.pdf

Links to Standards Conformance Statements for the Implementation

HL7

https://insiteone.com/s/UCA_R11_1_HL7_Conformance_Statement.pdf

DICOM

https://insiteone.com/s/UCA_R11_1_DICOM_Conformance_Statement.pdf

Link to general information on IHE: <https://www.ihe.net/>

IHE Integration Statement			
Vendor	Product Name	Version	Date
InsiteOne LLC	Unified Clinical Archive	11.1	8/20/2022

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profile	Actor	Integration Profile Option
Mammography Image	Image Manager/Archive	None
NM Image	Image Manager/Archive	None
Patient Demographics Query	Patient Demographics Consumer	None
Patient Information Reconciliation	Image Manager/Archive	None
Patient Information Reconciliation	Performed Procedure Step Manager	None
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
Simple Image and Numeric Report	Report Repository	None
Scheduled Workflow	Image Manager/Archive	None
Scheduled Workflow	Performed Procedure Step Manager	None
Scheduled Workflow.b	Image Manager/Archive	None
Scheduled Workflow.b	Performed Procedure Step Manager	None
Web Based Image Access	Image Document Responder	None
Web Based Image Access	Image Document Source	None

Internet address for vendor's IHE information: https://insiteone.com/s/UCA_R11_1_IHE_Integration_Statement.pdf	
Links to Standards Conformance Statements for the Implementation	
HL7	https://insiteone.com/s/UCA_R11_1_HL7_Conformance_Statement.pdf
DICOM	https://insiteone.com/s/UCA_R11_1_DICOM_Conformance_Statement.pdf
Link to general information on IHE: https://www.ihe.net/	